EN
DESCRIPTION
<p>Renishaw LaserAbutments™ are custom-made devices intended for aiding prosthetic functional and aesthetic rehabilitation. The abutments are manufactured centrally by Renishaw to designs created by a dental laboratory using a Renishaw scanner and Dental CAD software*. Designs for LaserAbutments can be created using a range of third-party scanners and CAD software packages*. In such circumstances, the manufacture by Renishaw of such customer-supplied designs is in accordance with the applicable essential requirements of the Medical Devices Directive (93/42/EEC) and under a quality management system that complies with ISO 13485:2003.</p> <p>* See the Renishaw website for the full and up-to-date list and instructions for each.</p> <p>LaserAbutments are available as screw retained or cement retained and come supplied with a polished emergence profile. This is chemically polished and the supra-gingival surface is fully prepared as part of the manufacturing process.</p>
INFORMATION FOR THE TECHNICIAN
<p>General</p> <p>Avoid excessive porcelain masses.</p> <p>Any grit blasting should be done with 50 µm alumina grit at 5 bar (73 psi) using a pencil nozzle.</p> <p>The screw(s) supplied with the LaserAbutments should not be used for laboratory purposes. For laboratory use, purchase additional screws from Renishaw or the implant OEM.</p> <p>Design</p> <p>Maximum abutment height is 20 mm.</p> <p>Maximum angulation from the implant direction is 30°.</p> <p>Finishing</p> <p>Emergence profile should be polished before returning to the clinician.</p> <p>Supra-gingival surface requires grit blasting before porcelain lay-up.</p> <p>Oxide firing should be done for 5 minutes at 950-980 °C with vacuum.</p> <p>Porcelain application (screw retained only)</p> <p>Maximum thickness of 1.5 mm.</p> <p>Use ISO 9693 compliant porcelains with firing temperatures up to 980 °C.</p> <p>The CoCr has a CTE (25 – 500 °C) of (14.0 – 14.5) × 10⁻⁶ K⁻¹, use porcelains with a similar CTE.</p> <p>Use according to manufacturer's instructions and recommendations.</p> <p>Apply opaque material in 2 firings – apply a thin layer followed by a second opaque layer.</p> <p>Wash part under running water before applying the next ceramic coating.</p> <p>Remove the porcelain mechanically only. Acid removal will cause corrosion of the metal.</p>
INDICATIONS
<ul style="list-style-type: none">Cement screw retained crown. Cement retained bridge. Suitable for anterior and posterior restorations. Intended for a range of dental implant systems.
CONTRAINDICATIONS
<ul style="list-style-type: none">Allergies to any of the materials used including the following: Co, Cr, Mo, W, Si, Fe, Mn. Implant and abutment offset angulation exceeding 30°. Patients suffering from bruxism or significant malocclusion.
⚠ WARNINGS AND PRECAUTIONS
<p>ⓘ The manufacturer's instructions for the implant and any tooling must be closely followed, alongside the information provided in these instructions for use.</p> <p>Ⓢ Renishaw LaserAbutments are customised specifically to one patient and are suitable for single use only.</p> <p>ⓘ The manufacturer's instructions for the screws provided by Renishaw with the LaserAbutments must be closely followed, alongside the information provided in these instructions for use (http://www.dess-abutments.com/pdf/documento-informativo-en.pdf).</p> <p>These instructions for use assume that both dental clinicians and laboratory technicians have been sufficiently trained in their respective professions to carry out the design and placement of the LaserAbutments.</p> <p>Clinicians and technicians should exercise their professional judgment to determine whether LaserAbutments are suitable for a particular patient. Any reshaping of the abutment is therefore a clinical determination.</p> <p>Clinicians and technicians are encouraged, as part of their professional development obligations, to regularly review Renishaw's latest developments in its range of dental products and their applications. Such information is available on our website or by talking to your Renishaw representative.</p> <p>Ensure close cooperation is maintained between the clinician and technician throughout the implant and abutment fitting process.</p> <p>Ensure threads and internal surfaces are clean and free of any foreign matter before sterilisation.</p> <p>⚠ The abutments are packed and delivered in a non-sterile state. The packaging is unsuitable for sterilisation purposes. They should be sterilised using a vacuum steriliser (EN 13060 type B). They should be fitted immediately following this and aseptic handling is imperative.</p> <p>The abutments should be used immediately after sterilisation.</p> <p>The usual clinical measures should be taken to avoid any swallowing or aspiration of the prosthetic components by the patient.</p> <p>Ensure that the patient's implant is fully integrated and healed before placement.</p> <p>Do not fit a pre-used screw for final placement of the abutment. Fit an unused screw – only Renishaw supplied screws to be used with this abutment.</p>
UNIVERSAL INFORMATION
<p>The abutments are supplied with screws that are specific to the abutment design submitted and can be secured using the implant OEM's driver. Follow the manufacturer's instructions and do not over-tighten or under-tighten the screws.</p>

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DE
INFORMATION FOR THE TECHNICIAN
<p>General</p> <p>Avoid excessive porcelain masses.</p> <p>Any grit blasting should be done with 50 µm alumina grit at 5 bar (73 psi) using a pencil nozzle.</p> <p>The screw(s) supplied with the LaserAbutments should not be used for laboratory purposes. For laboratory use, purchase additional screws from Renishaw or the implant OEM.</p> <p>Design</p> <p>Maximum abutment height is 20 mm.</p> <p>Maximum angulation from the implant direction is 30°.</p> <p>Finishing</p> <p>Emergence profile should be polished before returning to the clinician.</p> <p>Supra-gingival surface requires grit blasting before porcelain lay-up.</p> <p>Oxide firing should be done for 5 minutes at 950-980 °C with vacuum.</p> <p>Porcelain application (screw retained only)</p> <p>Maximum thickness of 1.5 mm.</p> <p>Use ISO 9693 compliant porcelains with firing temperatures up to 980 °C.</p> <p>The CoCr has a CTE (25 – 500 °C) of (14.0 – 14.5) × 10⁻⁶ K⁻¹, use porcelains with a similar CTE.</p> <p>Use according to manufacturer's instructions and recommendations.</p> <p>Apply opaque material in 2 firings – apply a thin layer followed by a second opaque layer.</p> <p>Wash part under running water before applying the next ceramic coating.</p> <p>Remove the porcelain mechanically only. Acid removal will cause corrosion of the metal.</p>
INDICATIONS
<ul style="list-style-type: none">Crown sealed or cemented. Bridge sealed. Convient aux restaurations antérieures et postérieures. Destiné à être utilisé avec une gamme de systèmes d'implants.
CONTRAINDICATIONS
<ul style="list-style-type: none">Allergy to any of the materials used, including the following: Co, Cr, Mo, W, Si, Fe, Mn. <p>In very extreme cases, electrochemical reactions or allergies to the contents of Renishaw CoCr DG1 may occur.</p> <ul style="list-style-type: none">Implant and abutment offset angulation exceeding 30°. Patients suffering from bruxism or significant malocclusion.
⚠ WARNINGS AND PRECAUTIONS
<p>ⓘ The manufacturer's instructions for the implant and any tooling must be closely followed, alongside the information provided in these instructions for use.</p> <p>Ⓢ Renishaw LaserAbutments are customised specifically to one patient and are suitable for single use only.</p> <p>ⓘ The manufacturer's instructions for the screws provided by Renishaw with the LaserAbutments must be closely followed, alongside the information provided in these instructions for use (http://www.dess-abutments.com/pdf/documento-informativo-en.pdf).</p> <p>These instructions for use assume that both dental clinicians and laboratory technicians have been sufficiently trained in their respective professions to carry out the design and placement of the LaserAbutments.</p> <p>Clinicians and technicians should exercise their professional judgment to determine whether LaserAbutments are suitable for a particular patient. Any reshaping of the abutment is therefore a clinical determination.</p> <p>Clinicians and technicians are encouraged, as part of their professional development obligations, to regularly review Renishaw's latest developments in its range of dental products and their applications. Such information is available on our website or by talking to your Renishaw representative.</p> <p>Ensure close cooperation is maintained between the clinician and technician throughout the implant and abutment fitting process.</p> <p>Ensure threads and internal surfaces are clean and free of any foreign matter before sterilisation.</p> <p>⚠ The abutments are packed and delivered in a non-sterile state. The packaging is unsuitable for sterilisation purposes. They should be sterilised using a vacuum steriliser (EN 13060 type B). They should be fitted immediately following this and aseptic handling is imperative.</p> <p>The abutments should be used immediately after sterilisation.</p> <p>The usual clinical measures should be taken to avoid any swallowing or aspiration of the prosthetic components by the patient.</p> <p>Ensure that the patient's implant is fully integrated and healed before placement.</p> <p>Do not fit a pre-used screw for final placement of the abutment. Fit an unused screw – only Renishaw supplied screws to be used with this abutment.</p>
UNIVERSAL INFORMATION
<p>The abutments are supplied with screws that are specific to the abutment design submitted and can be secured using the implant OEM's driver. Follow the manufacturer's instructions and do not over-tighten or under-tighten the screws.</p>

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